



Media for equine embryo transfer

(+ Designed to assure highest performance

- + Manufactured under GMP conditions in dedicated clean room
- + Gas-tight packaging for long shelf life
- + Convenient package sizes
- + Quality control of each batch performed under veterinarian supervision
- + Batch certificate available

Recovery - EquiFlush

Ready-to-use medium for equine embryo recovery. Contains buffers, BSA or PVA and antibiotics.

- Available in plastic bottles and infusion bags
- Bottles with rubber stopper for easy connection of the Y-tubing and tab on the bottom for hanging up
- Plastic bottles are made from only one material and are therefore recyclable
- Infusion bags with 3-way port configuration and sealed ports for hygienic fluid removal
- Graduation in steps of 100 ml simplifies dosage and allows monitoring the volume during flushing
- All media are ready to use (no additives required)

EquiFlush with BSA	
bag, 1 l	19982/6201
bag, 2 l	19982/6202
bottle, 1 l	19982/6205

EquiFlush with PVA	
bag, 1 l	19982/6211
bag, 2 l	19982/6212
bottle, 1 l	19982/6213



minitube



■ Holding - EquiHold

Ready-to-use holding medium for equine embryos. Contains buffers, glucose, pyruvate, BSA, amino acids, growth factors and antibiotics. Applicable for:

- Holding embryos at ambient conditions
- Transporting embryos
- Washing and rinsing embryos
- Loading embryos for transfer to recipient animal



EquiHold, 3 x 5 ml

19982/6250

■ Vitrification - Equine VIT-Kit

Vitrification is a proven method for cryopreservation of equine embryos. It causes less damage to the embryo than slow freezing and does not require a special freezing device. This vitrification system allows direct transfer after thawing of the embryo. No dilution steps are necessary after thawing.

The Equine VIT-Kit contains everything needed for vitrification of three embryos:

- Equilibration medium (3 x VS1, 3 x VS2, 3 x VS3)
- 3 x diluent
- 1 x syringe 1 ml for loading straws
- 1 x pipette tip (straw adapter)
- 4 x 5-well culture dish
- 10 x 0.25 ml sterilized straws for vitrification
- 10 x 0.5 ml sterilized straws for transfer
- 3 x ID rod
- 1 x aluminum cane
- 3 x goblet 10 mm



Equine VIT-Kit

19982/6251